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## INSTRUCTIONS FOR USE OF APG™ SUTURE

SYNTHETIC ABSORBABLE SURGICAL SUTURE

### DESCRIPTION

The APG™ Suture is a synthetic, absorbable, sterile, multifilament surgical suture, made of polyglycolic acid covered polycaprolactone and calcium stearate the APG™ Suture is available dyed violet to enhance its visibility on the tissue or on its natural color (beige). The pigment used in the violet suture is D&C #2 approved by the FDA. This product meets specifications set forth for Absorbable Surgical Sutures by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia (EP) and for the ABNT norm NBR 13904 - Filaments for Surgical Sutures. APG™ Suture is sterilized by ethylene oxide.

### INDICATION

The APG™ Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not indicated for use in cardiovascular or neurological tissues.

### ACTIONS

The APG™ Suture elicits a minimal acute inflammatory reaction in tissue followed by a gradual encapsulation of the suture by the fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of suture APG™ occurs by means of hydrolysis, where the polymer degrades to glycolic acid which is subsequently absorbed and metabolized in the body. Absorption begins as loss of tensile strength followed by a loss of mass. Absorption studies in vitro indicate that after 14 days the tensile strength is approximately 60% of its initial tensile strength. The suture absorption is completed between 90-120 days after the implant.

### CONTRAINDICATION

The APG™ Suture being absorbable should not be used where extended approximation of tissues under stress is required.

This device is contraindicated to persons that could have any allergic reaction to the medical device material.

### WARNINGS

Users must be acquainted with surgical procedures and practices involving absorbable sutures before using APG™ Suture to close incisions for the risk of wound dehiscence may vary according to the site of application and the suture material used. The use of this suture may be unsuitable for elder, undernourished, weak or other patients whose conditions may slow down the wound from skinning over. As this is an absorbable suture material the use of supplemental non absorbable sutures should be considered by the surgeon in the closure of wounds which may undergo expansion, stretching or distension, or which may require additional support.

The APG™ Suture as an absorbable suture it may act transiently as a foreign body. As with any foreign body prolonged contact of any suture with saline solution, such as those found at the biliary and urinary canals, may result in calculi formation. The suture should only be used for or under surgeon's orientation according to applicable surgical techniques and it should only be manipulated by qualified people. Do not use the suture for other purposes, just the one indicated by the manufacturer. When handling this suture be sure of doing it aseptically. Acceptable surgical practices must be followed to treat any infected or contaminated wounds. Do not re-sterilize. Discard the open packages and unused sutures. Discard left over. Do not use if packing is violated or damp. Do not use the product beyond the validity period. Disposal of used sutures: Risk of biological contamination - adopt internal disposal procedures, taking into account the risk of biological contamination. Disposal of used and non-used sutures: Perforating and cutting materials - adopt internal disposal procedures for perforating and cutting materials.

### PRECAUTIONS

Under certain circumstances, especially in orthopedic procedures, the immobilization of joints by external support may be employed, if the surgeon finds it necessary. It is necessary to reflect on the use of absorbable sutures in tissues with low supply of blood, for it may occur extrusion of the suture and the slowing down of the absorption. In the handling this or any other suture material care should be taken to avoid damage from handling. Avoid the crimping and the crushing of the needles by the inadequate handling of tweezers and needle holders. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens. Dispose used needles in recipients for potentially contaminated material.

### USE

Use the necessary amount of suture according to the wound extension and the technique adopted by the surgeon. The size and diameter of suture and needle used for each type of surgical intervention are defined as per medical criteria.

### ADVERSE REACTIONS

Adverse reactions connected with the use of this absorbable synthetic suture include: wound dehiscence, infection, tissue minimal acute inflammatory reaction, localized irritation, suture extrusion and slow absorption in tissues with poor supply of blood. Calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs. The surgeon must inform the patient to look for medical assistance in case of any suspect of adverse reactions symptoms.

### STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30°C (59°F to 86°F), protected from direct sunlight, heat sources, humidity and dust.

### HOW SUPPLIED

The APG™ Suture is supplied dyed in violet or in its natural color (beige), in sterile envelopes, in U.S.P. sizes 8-0 through (metric sizes 0.4 through 5). Suture is supplied in a variety of lengths, attached to many types of non traumatic needles, according to each type of surgical intervention, in boxes with 12, 24 or 36 units.

### EXPIRATION

Five years from manufacturing date. See box and envelope.

PACKAGES WITH 12 (TWELVE), 24 ( TWENTY FOUR) OR

36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.

SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES

DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN

DESTROY AFTER USING

STERILIZATION METHOD : ETHYLENE OXIDE

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CNPJ 02.370.649/0001-20 - Fabrica: Tel: 55 19 3646-7300

Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

ANVISA REGISTRATION N.º POLYGLYCOLIC ACID WITH NEEDLE: 10372360004